

## **Installation**

1. Receive system components at dentist's office.
2. Check to ensure network cable has been run. If not, do not proceed until it has.
3. Once network cable has been run, place each system component in the room it will be connected.
4. Unpack each system component and place in proper location.
5. Install power to each system component.
6. Connect monitors, keyboards, and mouses to corresponding computer.
7. Connect network printer and each computer to corresponding wall plate network connection.
8. Connect system component LAN cable feeds from wall plate to router.
9. Connect hi-speed modem to phone line and router.
10. Installation is complete.

## **Testing**

1. Power up each system component.
2. Start up practice management software on each computer to verify it has been installed.

3. Start up digital radiography software on each operatory computer to verify it has been installed.
4. From each front desk computer, search for several patients in practice management software to ensure patient database is properly installed and connected to software.
5. From each operatory computer, call up patient in practice management software and click on link to digital radiography software to ensure functionality.
6. Hook up x-ray sensor to each operatory computer and capture test x-ray to ensure functionality.
7. Testing is complete.